

# **PowerRef 3**

#### SPECIFICATIONS



Plusoptix GmbH Neumeyerstr. 48 90411 Nuremberg Germany

 Phone:
 +49-911-59 83 99-10

 Fax:
 +49-911-59 83 99-90

 E-Mail:
 c.schmidt@plusoptix.com

 Internet:
 www.plusoptix.com



# **Technical specification of PowerRef 3**

#### Measurements

Refraction (vertical meridian)	binocular and monocular	
Spherical range	+5.00/-7.00 dpt in 0.01 dpt steps	
Pupil size	4.0 – 8.0 mm in 0.1 mm steps ± 10%	
Pupil distance	mm in 1.0 mm steps ± 10 %	
Time per measurement	0.02 s (i.e. 50 Hz sampling rate)	
Measuring distance	1 m ( ± 5 cm)	
Measurement principal	binocular, dynamic photosciascopy	

#### Interfaces and standards

Interfaces 4 x USB, 1 x VGA (1024 x 769)	
--	--

#### Power

Main 100 - 240 VAC (50 -	- 60Hz), 2.5 A (excl. monitor)
--------------------------	--------------------------------

### Environmental requirements

Temperature	operating	10 to 35 ℃ (50 to 92 ℉)
	storage	0 to 50 ℃ (32 to 122 ℉)
Humidity	operating	20 - 80 % not condensing
	storage	10 - 85 % not condensing

## **Documentation specification of PowerRef 3**

PowerRef 3 readings are saved in a csv-file. This file is called "logfile" and contains the following information:

Value	Description
Time [ms]	Time stamp of the measured values. Every 20ms there will be one measured value.
Brightness [%]	LED brightness in percent
PupilFound (L and R)	1: Pupil found; 0: Pupil not found.
PupilsizeMM(L and R).X() [mm]	Pupil width right eye and left eye
PupilsizeMM(L and R).Y() [mm]	Pupil height right eye and left eye
PupilBrightness(L and R) [DU]	Qualitative pupil brightness left eye and right eye in values between 0 and 255.
PurkinjeSize(L and R).X() [px]	Horizontal dimension of the purkinje reflex in pixel.
PurkinjeSize(L and R).Y() [px]	Vertical dimension of the purkinje reflex in pixel.
Decentration(L and R).X() [px]	Horizontal distance of the purkinje reflex from the middle of the pupil.
Decentration(L and R).Y() [px]	Vertical distance of the purkinje reflex from the middle of the pupil.
Gaze(L and R).X() [¶	Angle of the horizontal gaze in degrees.
Gaze(L and R).Y() [¶	Angle of the vertical gaze in degrees.
Refraction(L and R).Sphere() [dpt]	Measured refraction in diopters of each frame.
RefractionLeft(resp. Right) Contents()	Value between 0-10. Describes the actual state of the average calculation of the refractionLeft(resp. Right).Mean() [dpt].
RefractionLeft(resp. Right).Mean() [dpt]	Average of the last 10 frames in diopters.
RefractionLeft(resp. Right).StdDeviation() [dpt]	The standard deviation of the last ten frames.
InterpupDist() [mm]	Pupil distance in mm.
InterpupAngle() [ <sup>9</sup>	Angle of the camera to the hori zontal.



## **PowerRef 3 set-up schematics**

Using a hot mirror, PowerRef 3 can be moved out of patient's line of sight. This way visual targets can be presented at any distance.

